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# User Manual for Cleaning and Sterilization of Dental Implant Instruments

## User Manual for Cleaning Dental Implant Instruments

Cowellmedi Co., Ltd. recommends using an enzymatic detergent solution, a user-friendly product, generating a bit of foams and containing a surfactant and easily cleansing pollutants.

Only the detergent solution having the proved effects (FDA-approved, CE mark) in selecting an enzymatic detergent solution should be used.

Because various detergent solutions and disinfectants exist in the world, Cowellmedi does not recommend any specific brand.

The detergent solution used in our user manual for cleaning dental implant instruments is the one whose effectiveness has been proved. The detergent solution is as follows:

CIDEZYME (ENZYMATIC DETERGENT SOLUTION): J&J

All detergent solutions can generate their optimal performance effects, only if the detergent solution's dilution ratio and cleaning temperature recommended by the manufacturer concerned are prepared.

**\*Note** Since cleaning variables may occur according to detergent equipment or human factors in manual cleaning carried out in medical facilities, proper procedure or document, whose effectiveness has been proved, should be in place.

## User Manual for Sterilization of Dental Implant Instruments

The sterilization guideline of dental implant instruments recommended by Cowellmedi Co., Ltd. is the high pressure steam sterilization method, whose effectiveness has been proved, and 10-6 sterilization assurance level (SAL) is assured.

A sterilizer used in a medical facility needs to comply with the recommendation of the manufacturer. Upon sterilizing several instruments, the specifications of the sterilizer recommended by the manufacturer concerned should not be exceeded.

An implant instrument set should be put on a tray, made for it to directly contact the instrument surface, with easy steam penetration and in a separate case for sterilization preparation.

**\*Note** Because sterilization variables may occur according to sterilization equipment or human factors in terms of sterilization handled in medical facilities, proper procedure or document, whose effectiveness has been proved, should be in place.

### ►How to clean an implant instrument

The guideline to clean an implant instrument recommended by Cowellmedi presents the cleaning method, whose effectiveness has been proved, and the cleaning method is presented below:

1. Soak an instrument in sufficient amount of ethyl alcohol (70%), and carry out ultrasonic cleaning.
2. To remove bloodstains and alien substances that can remain, even after the cleaning with ethyl alcohol, clean the parts, difficult to be cleaned with distilled water or flowing water, using a syringe or a pipe cleaner (soft brush for cleaning).
3. After soaking the instrument in sufficient water to fully sock it by diluting the enzymatic detergent solution with the tap water or distilled water in the ratio recommended by the manufacturer, carry out ultrasonic cleaning at 35~40kHz for 15 minutes.
4. To remove the enzymatic detergent solution used to clean the instrument and alien substances, rinse with flowing distilled water for 3 minutes.
5. After completely removing moisture with dry cloth or hot air blower, perform sterilization in a user sterilization method by wrapping the assembled instrument on a tray or the instrument alone with surgical drape.

### ►Cautions

1. Be careful not to lose the parts of the dental implant instrument in cleaning the instrument.
2. Be careful not be hurt by sharp part of the dental implant instrument in handling the instrument.
3. We guarantee a number of reusable number for the COWELL® Reusable Surgical Instrument within 60 cycles.

### ►How to sterilize dental implant instruments

The user manual for sterilization of dental implant instruments recommended by Cowellmedi presents a sterilization method, whose effectiveness has been proved, and the method is as follows;

The sterilization method and conditions is moist heat sterilization method (high pressure steam sterilization). Cowellmedi recommends sterilization at **121℃ for 40min or more** by putting the instrument on its tray and wrapping it with surgical drape.

### ►Cautions

1. Check whether any abnormality such as alien substances and rust exist in the instrument before using the instrument.
2. The instrument is supplied in the non-sterilized state, and therefore sterilization needs to be carried out before using it.
3. We recommend immediately cleaning the instrument according to manufacturer's cleaning method, after it is used. After cleaning, the instrument needs to be stored in a dry and clean place. Before using the instrument, sterilization must be carried out.
4. We guarantee a number of reusable number for the COWELL® Reusable Surgical Instrument within 60 cycles.

### Cowellmedi Co., Ltd.

48 Hakgam-daero, 221 beon-gil, Sasang-gu, Busan 46986, Republic of Korea Tel. +82-51-312-2027-8 Fax. +82-51-316-2628 <http://www.cowellmedi.com>