COWELL REGENERATIVE SOLUTION

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
- 1 Analysis of hydrolyzable polyethylene glycol hydrogels and deproteinized bone mineral as delivery systems for glycosylated and non-glycosylated bone morphogenetic protein-2. *Acta Biomater.* 2012 Jan;8(1):116-23.
- 2 Effects of rhBMP-2 Coating Tricalcium Phosphate on Socket Preservation in Dog Extraction Socket. *Tissue Engineering and Regenerative Medicine*, Vol. 5, No. 4-6, pp 637-642 (2008)
- 3 Effects of Polycaprolactone-Tricalcium Phosphate, Recombinant Human Bone Morphogenetic Protein-2 and Dog Mesenchymal Stem Cells on Bone Formation: Pilot Study in Dogs. *Yonsei Med J* 50(6): 825-831, (2009).
- 4 The induction of bone formation in rat calvarial defects and subcutaneous tissues by recombinant human BMP-2, produced in *Escherichia coli*. *Biomaterials* 31 (2010) 3512-3519.
- 5 Alveolar ridge augmentation using anodized implants coated with *Escherichia coli*-derived recombinant human bone morphogenetic protein 2. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod.* (2011) Jul;112(1):42-9.
- 6 Bone formation of *Escherichia coli* expressed rhBMP-2 on absorbable collagen block in rat calvarial defects. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2011;111:298-305.
- 7 Bone formation of block and particulated biphasic calcium phosphate lyophilized with *Escherichia coli*-derived recombinant human bone morphogenetic protein 2 in rat calvarial defects. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2011;112:298-306.
- 8 Induction of bone formation by *Escherichia coli*-expressed recombinant human bone morphogenetic protein-2 using block-type macroporous biphasic calcium phosphate in orthotopic and ectopic rat models. *J Periodontal Res.* (2011) Dec; 46(6):682-90.
- 9 Enhanced adipogenic differentiation and reduced collagen synthesis induced by human periodontal ligament stem cells might underlie the negative effect of recombinant human bone morphogenetic protein-2 on periodontal regeneration. *J Periodontal Res* (2011); 46: 193-203.
- 10 The Effects of rhBMP-2 Injection at Distraction Osteogenesis of Rats' Tibia. *Tissue Engineering and Regenerative Medicine*, Vol. 8, No. 2, pp 158-163 (2011).
- 11 Discontinuous Release of Bone Morphogenetic Protein-2 Loaded Within Interconnected Pores of Honeycomb-Like Polycaprolactone Scaffold Promotes Bone Healing in a Large Bone Defect of Rabbit Ulna. *Tissue Eng Part A.* 2011 Oct;17(19-20):2389-97.
- 12 The effect of immobilization of heparin and bone morphogenetic protein-2 to bovine bone substitute on osteoblast-like cell's function. *J Adv Prosthodont* 2011; 3:145-51.
- 13 Multicenter, randomized clinical trial on the efficacy and safety of *Escherichia coli*-derived rhBMP-2 with β -Tricalcium phosphate and hydroxyapatite in human extraction sockets. *J Adv Prosthodont* 2011; 4:178-182.
- 14 Effects of Anodized Implants Coated With *Escherichia coli*-Derived Recombinant Human Bone Morphogenetic Protein-2 on Osseointegration in Rabbits. *Tissue Engineering and Regenerative Medicine*, Vol. 8, No. 1, pp 62-68 (2011).
- 15 Novel analysis model for implant osseointegration using ectopic bone formation via the recombinant human bone morphogenetic protein-2/macroporous biphasic calcium phosphate block system in rats: a proof-of concept study. *J Periodontal Implant Sci* 2012; 42:136-143.
- 16 Effects of anodized implants coated with *Escherichia coli*-derived rhBMP-2 in beagle dogs. *Int. J. Oral Maxillofac. Surg.* 2012; 41: 1577-1584.
- 17 Bone formation of middle ear cavity using biphasic calcium phosphate lyophilized with *Escherichia coli*-derived recombinant human bone morphogenetic protein 2 using animal model. *International Journal of Pediatric Otorhinolaryngology* 77 (2013) 1430-1433.
- 18 Bone formation and remodeling of three different dental implant surfaces with *Escherichia coli*-derived recombinant human bone morphogenetic protein 2 in a rabbit model. *Int J Oral Maxillofac Implants.* 2013; 28(2):424-30.
- 19 Recombinant Human Bone Morphogenetic Protein-2 Stimulates the Osteogenic Potential of the Schneiderian Membrane: A Histometric Analysis in Rabbits. *Tissue Eng Part A.* 2013 Sep;19(17-18):1994-2004.
- 20 The effect of anodized implants coated with combined rhBMP-2 and recombinant human vascular endothelial growth factors on vertical bone regeneration in the marginal portion of the peri-implant. *Oral Surg Oral Med Oral Pathol Oral Radiol* 2013;115:e24-e31.
- 21 Sinus augmentation using BMP-2 in a bovine hydroxyapatite/collagen carrier in dogs. *J Clin Periodontol* 2014; 41: 86-93.
- 22 Low-Dose Recombinant Human Bone Morphogenetic Protein-2 to Enhance the Osteogenic Potential of the Schneiderian Membrane in the Early Healing Phase: In Vitro and In Vivo Studies. *J Oral Maxillofac Surg* 72:1480-1494, 2014.
- 23 Prospective randomized, controlled trial of sinus grafting using *Escherichia coli*-produced rhBMP-2 with a biphasic calcium phosphate carrier compared to deproteinized bovine bone. *Clin Oral Implants Res.* 2015 Dec;26(12):1361-8.
- 24 Controlled release of BMP-2 using a heparin-conjugated carrier system reduces in vivo adipose tissue formation. *J Biomed Mater Res A.* 2015 Feb;103(2):545-54.
- 25 The efficacy of BMP-2 preloaded on bone substitute or hydrogel for bone regeneration at peri-implant defects in dogs. *Clin Oral Implants Res.* 2015 Dec;26(12):1456-65.
- 26 Effect of rhBMP-2 Immobilized Anorganic Bovine Bone Matrix on Bone Regeneration. *Int. J. Mol. Sci.* 2015, 16, 16034-16052.
- 27 Effects of rhBMP-2 on Sandblasted and Acid Etched Titanium Implant Surfaces on Bone Regeneration and Osseointegration: Split-Mouth Designed Pilot Study. *Biomed Res Int.* 2015; 2015:459393.
- 28 Comparison of collagen membrane and bone substitute as a carrier for rhBMP-2 in lateral onlay graft. *Clin Oral Implants Res.* 2015;26(1):e13-9.
- 29 Effects of BMP-2 Delivery in Calcium Phosphate Bone Graft Materials with Different Compositions on Bone Regeneration. *Materials* 2016, 9, 954.
- 30 Source and Carrier Effect on the Bioactivity of BMP Bio-Implants. *Master of Science* 2013. Sylvie Di Lullo 2013, Faculty of Dentistry, University of Toronto.
- 31 Soft and hard tissue changes when socket preservation using rhBMP-2, PRP and Non-Resorbable dPTFE membrane. *Dental Implant Journal: Vol. 3, May, 2014.*
- 32 The effect of rhBMP-2 bonegraft on infrabony defects. *Dental Implant Journal: Vol. 3, May, 2014.*

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Clinical Case I


Fig 01. Preoperative radiograph.


Fig 02. Preoperative CBCT image
Sinusitis in both sinus cavities.

Fig 03. Incision and flap elevation.
Removal of granulation tissue.

Fig 04. Suction of pus from the sinus.


Fig 05. Bone grafting with InnoGraft B.
Resorbable membrane application.


Fig 06. Postoperative radiograph.

Fig 07. Postoperative CBCT image.


Fig 08. Postoperative radiograph at
6 weeks.
Final restoration delivery.


Fig 09. CBCT image of postoperative 10
months.

Clinical Case II


Fig 01. Preoperative radiograph. 3 months after extraction in lower left posterior.


Fig 02. Incision and flap elevation.


Fig 03. Implant placement on #35, 37.

Fig 04. Bone grafting with InnoGraft B and InnoOss Allo.

Fig 05. Non-resorbable membrane application.

Fig 06. Sture.


Fig 07. Postoperative radiograph.


Fig 08. Postoperative CBCT image of #35 (Lt), #37 (Rt).

Fig 09. Clinical view of postoperative 2 weeks.

Fig 10. Clinical view of postoperative 4 months. Final restoration delivery.

Fig 11. Radiograph of postoperative 4 months.

Fig 15. 4 months postoperative CBCT image of #35 (Lt), #37 (Rt).

Clinical Case III


Fig 01. Preoperative radiograph.

Fig 02. Preoperative clinical view.


Fig 03. Extraction and immediate implant placement on #14, #15.


Fig 04. Bone grafting with InnoGraft B.


Fig 05. Suture.


Fig 06. Postoperative radiograph.


Fig 07. Postoperative CBCT image of #14 (Lt), #15 (Rt).

Fig 08. Clinical view of postoperative 3 months. Abutment connection.

Fig 09. Final restoration delivery.


Fig 10. Radiograph of postoperative 3 months.


Fig 11. 3 months postoperative CBCT image of #14 (Lt), #15 (Rt).

Clinical Case IV

Fig 01. Preoperative radiograph.
Horizontal fracture of #13.

Fig 02. Preoperative clinical view.
Old prosthesis removal.


Fig 03. Extraction and immediate implant placement on #13, #14.


Fig 04. Bone grafting with InnoGraft B.


Fig 05. Postoperative radiograph.

Fig 06. Postoperative CBCT image #13.

Fig 07. Clinical view of postoperative 4 months. Abutment connection.

Fig 08. Final restoration delivery.


Fig 09. Radiograph of postoperative 4 months.

Fig 10. 4 months postoperative CBCT image of #13.

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